

DATA SHEET OVERVIEW I

SPIROVITAL THERAPY/GASTROVITAL THERAPY


LITTLE ATMOS



LITTLE ATMOS PLUS



BASIS PLUS



The little power dwarf

- 1 Original A+ Bicat
- Integrated colour light therapy (blue)
- 3 intensity levels (50 %, 75 %, 100 %)
- Microcontroller-controlled electronics
- LCD Display
- Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES)
- Individual programme settings
- Automatic self-test
- Switchable acoustic signal
- Main device switch & integrated port (DC) for Stream & Stream HS incl. switch

Optionally available: Aroma therapy

The power dwarf with a big plus

- 1 Original A+ Bicat, 1 A+ Single/Classic-Catalyst
- Integrated colour RGB light therapy (8 colours)
- 3 intensity levels (50 %, 75 %, 100 %)
- Microcontroller-controlled electronics
- LCD Display
- Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES)
- Individual programme settings
- Automatic self-test
- Switchable acoustic signal
- Main device switch & integrated port (DC) for Stream & Stream HS incl. switch

Aroma therapy
(included in the scope of delivery)

The perfect start or as a second device

- 2 Original A+ Single/Classic-Catalysts
- Integrated colour light therapy (blue)
- 3 intensity levels (50 %, 75 %, 100 %)
- Microcontroller-controlled electronics
- LCD Display
- Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES)
- Individual programme settings
- Automatic self-test
- Switchable acoustic signal
- Integrated port (DC) for Stream & Stream HS incl. switch

Optionally available: Aroma therapy

Power consumption	max. 1,185 mAh
Power supply	100-240V AC 50-60Hz Out: 12V DC 2.5 A
Dimensions (mm) incl. glass bottle	210 x 235 x 201 (W x H x D)
Weight incl. accessories	2.4 kg

Power consumption	max. 1,440 mAh
Power supply	100-240V AC 50-60Hz Out: 12V DC 2.5 A
Dimensions (mm) incl. glass bottle	210 x 235 x 201 (W x H x D)
Weight incl. accessories	2.7 kg

Power consumption	max. 1,640 mAh
Power supply	100-240V AC 50-60Hz Out: 12V DC 2.5 A
Dimensions (mm) incl. glass bottle	345 x 305 x 240 (W x H x D)
Weight incl. accessories	3.8 kg

Service	Every 2,000 hours of operation (approx. 6,000 sessions)
Warranty	2 years
Shipment includes	Power supply unit, Nasal cannulas (4x), Nasal cannulas for children (1x), Glass bottle incl. sparkling element (1x), Air filters V2 (2x), Nasal cannula bag, User manual, Quick guide

Article-No. (EAN-No.) [2022055 \(079297116773\)](#)

Service	Every 2,000 hours of operation (approx. 6,000 sessions)
Warranty	2 years
Shipment includes	Power supply unit, Nasal cannulas (4x), Nasal cannulas for children (1x), Glass bottle incl. sparkling element (1x), Air filters V2 (2x), Nasal cannula bag, User manual, Quick guide, Aroma bottle incl. oil (1x)

Article-No. (EAN-No.) [2021055 \(0792971115820\)](#)

Service	Every 2,000 hours of operation (ca. 3.000 Anwendungen)
Warranty	2 years
Shipment includes	Power supply unit, Nasal cannulas (10x), Glass bottle incl. sparkling element (1x), Air filters V2 (2x), Nasal cannula bag, User manual


Article-No. (EAN-No.) [2012055 \(0792971115813\)](#)



DATA SHEET OVERVIEW II

SPIROVITAL THERAPY/GASTROVITAL THERAPY

PROFESSIONAL PLUS compact



AvantGarde



TRAVEL PLUS V2



<p>The classic model for ambitious people 2 Original A+ Bicats</p> <p>Integrated colour RGB light therapy (8 colours)</p> <p>3 intensity levels (50 %, 75 %, 100 %) Microcontroller-controlled electronics LCD Display Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES) Individual programme settings Automatic self-test Switchable acoustic signal Integrated port (DC) for Stream & Stream HS incl. switch</p> <p>Optionally available: Aroma therapy</p>	<p>The design edition 2 Original A+ Bicats (platin, blanc, noir) 4 Original A+ Single/Classic-Catalysts (cristal) Integrated colour RGB light therapy (8 colours) Underbody lighting (nur cristal, blanc, platin) 3 intensity levels (50 %, 75 %, 100 %) Microcontroller-controlled electronics LCD Display Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES) Individual programme settings Automatic self-test Switchable acoustic signal</p> <p>Optionally available: Aroma therapy</p>	<p>For on the go 2 Original A+ Bicats</p> <p>Integrated colour RGB light therapy (8 colours)</p> <p>3 intensity levels (50 %, 75 %, 100 %) Microcontroller-controlled electronics LCD Display, status LED, charging indicator in display Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES) Individual programme settings Automatic self-test Switchable acoustic signal, TSA lock Device main switch & integrated port (DC) for Stream & Stream HS incl. switch</p> <p>Optionally available: Aroma therapy</p>
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<p>Power consumption max. 2,150 mAh</p> <p>Power supply 100-240V AC 50-60Hz Out: 12V DC 2.5 A</p> <p>Dimensions incl. glass bottle Professional Plus compact Professional Plus (classic)* 345 x 305 x 240 mm 345 x 305 x 340 mm</p> <p>Weight incl. accessories Professional Plus compact 4.1 kg Professional Plus (classic)* 4.7 kg</p>	<p>Power consumption cristal, platin, blanc noir* max. 3,150 mAh max. 2,150 mAh</p> <p>Power supply 100-240V AC 50-60Hz Out: 12V DC 3.8 A</p> <p>Dimensions (mm) incl. glass bottle 350 x 315 x 350 (W x H x D)</p> <p>Weight incl. accessories 7.95 kg</p>	<p>Power consumption max. <4,000</p> <p>Charger 100-240V AC 50-60Hz Out: 12,6V DC </=2.5A</p> <p>Device battery 12.6V / 7.75Ah</p> <p>Dimensions (mm) incl. glass bottle 350 x 315 x 350 (W x H x D)</p> <p>Weight incl. accessories 7.95 kg</p>
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<p>Service Every 2,000 hours of operation (approx. 6,000 sessions)</p> <p>Warranty 2 years</p> <p>Shipment includes Power supply unit, Nasal cannulas (20x), Glass bottle incl. sparkling element (1x), Air filters V2 (2x), Nasal cannula bag, User manual</p> <p>Article-No. (EAN-No.) Professional Plus compact 2040155 (0792971115837) Professional Plus (classic)* 2040055 (0792971115820)</p>	<p>Service Every 2,000 hours of operation (approx. 6,000 sessions)</p> <p>Warranty 2 years</p> <p>Shipment includes Power supply unit, Nasal cannulas (20x), Glass bottle incl. sparkling element (1x), Air filters V2 (2x), Nasal cannula bag, User manual</p> <p>Article-No. (EAN-No.) cristal 2080055 (0792971115844) platin 2081055 (0792971115851) blanc 2082055 (0792971115868) noir 2083055 (0792971115875)</p>	<p>Service Every 2,000 hours of operation (approx. 6,000 sessions)</p> <p>Warranty 2 years</p> <p>Shipment includes Charger, Battery (1x), Nasal cannulas (20x), Glass bottle incl. sparkling element (1x), Air filters V2 (2x), Nasal cannula bag, User manual, HS Diffusor (1x)</p> <p>Optional Up to 2 additional batteries for up to 1,000 min. operating time</p> <p>Article-No. 2021055</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



DATA SHEET OVERVIEW III

SPIROVITAL THERAPY/GASTROVITAL THERAPY

AvantGarde
DESIGN EDITION

AIRNERGY+
ENERGY FROM INSIDE

 <p>AvantGarde cristal</p>	 <p>AvantGarde blanc</p>	 <p>AvantGarde noir</p>	 <p>AvantGarde platin</p>
<p>4 Original A+ Single/Classic-Catalysts Integrated colour RGB light therapy (8 colours) Underbody lighting 3 intensity levels (50 %, 75 %, 100 %) Microcontroller-controlled electronics LCD Display Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES) Individual programme settings Automatic self-test Switchable acoustic signal</p> <p>Optionally available: Aroma therapy</p> 	<p>2 Original A+ Bicats Integrated colour RGB light therapy (8 colours) Underbody lighting 3 intensity levels (50 %, 75 %, 100 %) Microcontroller-controlled electronics LCD Display Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES) Individual programme settings Automatic self-test Switchable acoustic signal</p> <p>Optionally available: Aroma therapy</p> 	<p>2 Original A+ Bicats Integrated colour RGB light therapy (8 colours) 3 intensity levels (50 %, 75 %, 100 %) Microcontroller-controlled electronics LCD Display Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES) Individual programme settings Automatic self-test Switchable acoustic signal</p> <p>Optionally available: Aroma therapy</p> 	<p>2 Original A+ Bicats Integrated colour RGB light therapy (8 colours) Underbody lighting 3 intensity levels (50 %, 75 %, 100 %) Microcontroller-controlled electronics LCD Display Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES) Individual programme settings Automatic self-test Switchable acoustic signal</p> <p>Optionally available: Aroma therapy</p> 
<p>Power consumption max. 3,150 mAh</p> <p>Power supply 100-240V AC 50-60Hz</p> <p>Out: 12V DC 3.8A</p> <p>Dimensions (mm) incl. glass bottle 350 x 315 x 350 (W x H x D)</p> <p>Weight incl. 7.95 kg</p>	<p>Power consumption max. 3,150 mAh</p> <p>Power supply 100-240V AC 50-60Hz</p> <p>Out: 12V DC 3.8 A</p> <p>Dimensions (mm) incl. glass bottle 350 x 315 x 350 (W x H x D)</p> <p>Weight incl. 7.95 kg</p>	<p>Power consumption max. 2,150 mAh*</p> <p>Power supply 100-240V AC 50-60Hz</p> <p>Out: 12V DC 3.8 A</p> <p>Dimensions (mm) incl. glass bottle 350 x 315 x 350 (W x H x D)</p> <p>Weight incl. 7.95 kg</p>	<p>Power consumption max. 3,150 mAh</p> <p>Power supply 100-240V AC 50-60Hz</p> <p>Out: 12V DC 3.8 A</p> <p>Dimensions (mm) incl. glass bottle 350 x 315 x 350 (W x H x D)</p> <p>Weight incl. 7.95 kg</p>
<p>Service Every 2,000 hours of operation (approx. 6,000 sessions)</p> <p>Warranty 2 years</p> <p>Shipment includes Power supply unit, Nasal cannulas (20x), Glass bottle incl. sparking element (1x), Air filters V2 (2x), Nasal cannula bag, User manual</p> <p>Article-No. <u>2080055</u> (EAN-No.) <u>(0792971115844)</u></p>	<p>Service Every 2,000 hours of operation (approx. 6,000 sessions)</p> <p>Warranty 2 years</p> <p>Shipment includes Power supply unit, Nasal cannulas (20x), Glass bottle incl. sparking element (1x), Air filters V2 (2x), Nasal cannula bag, User manual</p> <p>Article-No. <u>2082055</u> (EAN-No.) <u>(0792971115868)</u></p>	<p>Service Every 2,000 hours of operation (approx. 6,000 sessions)</p> <p>Warranty 2 years</p> <p>Shipment includes Power supply unit, Nasal cannulas (20x), Glass bottle incl. sparking element (1x), Air filters V2 (2x), Nasal cannula bag, User manual</p> <p>Article-No. <u>2081055</u> (EAN-No.) <u>(0792971115851)</u></p>	<p>Service Every 2,000 hours of operation (approx. 6,000 sessions)</p> <p>Warranty 2 years</p> <p>Shipment includes Power supply unit, Nasal cannulas (20x), Glass bottle incl. sparking element (1x), Air filters V2 (2x), Nasal cannula bag, User manual</p> <p>Article-No. <u>2083055</u> (EAN-No.) <u>(0792971115875)</u></p>

*no underbody lighting

AIRNERGY_Data_sheet_Overview_2024



LITTLE ATMOS



LITTLE ATMOS PLUS



The little power dwarf

- 1 Original A+ Bicat
- Integrated colour light therapy (blue)
- 3 intensity levels (50 %, 75 %, 100 %)
- Microcontroller-controlled electronics
- LCD Display
- Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES)
- Individual programme settings
- Automatic self-test
- Switchable acoustic signal
- Main device switch & integrated port (DC) for Stream & Stream HS incl. switch
- Standard design (white); optionally available in 9 different designs

Optionally available: Aroma therapy
 Special or series production (OEM) possible on request



The power dwarf with a big plus

- 1 Original A+ Bicat, 1 Original A+ Single/Classic-Catalyst
- Integrated colour RGB light therapy (8 colours)
- 3 intensity levels (50 %, 75 %, 100 %)
- Microcontroller-controlled electronics
- LCD Display
- Multi-lingual menu guidance (DE, EN, NL, FR, IT, ES)
- Individual programme settings
- Automatic self-test
- Switchable acoustic signal
- Integrated port (DC) for Stream & Stream HS incl. switch
- Power switch
- Standard design (white); optionally available in 9 different designs

Aroma therapy (included in the scope of delivery)
 Special or series production (OEM) possible on request



Power consumption	max. 1,185 mAh
Power supply	100-240V AC 50-60Hz Out: 12V DC 2.5 A
Dimensions (mm) incl. glass bottle	210 x 235 x 201 (W x H x D)
Weight incl. accessories	2.4 kg

Power consumption	max. 1,440 mAh
Power supply	100-240V AC 50-60Hz Out: 12V DC 2.5 A
Dimensions (mm) incl. glass bottle	210 x 235 x 201 (W x H x D)
Weight incl. accessories	2.7 kg

Service	Every 2,000 hours of operation (approx. 6,000 sessions)
Warranty	2 years
Shipment includes	Power supply unit, Nasal cannulas (4x), Nasal cannulas for children (1x), Glass bottle incl. sparkling element (1x), Air filters V2 (2x), Nasal cannula bag, User manual, Quick guide
Article-No. (EAN-No.)	2022055 (079297116773)

Service	Every 2,000 hours of operation (approx. 6,000 sessions)
Warranty	2 years
Shipment includes	Power supply unit, Nasal cannulas (4x), Nasal cannulas for children (1x), Glass bottle incl. sparkling element (1x), Air filters V2 (2x), Nasal cannula bag, User manual, Aroma bottle incl. oil (1x), Quick guide
Article-No. (EAN-No.)	2021055 (0792971115820)



STREAM silver STREAM carbon



STREAM HS V2







Allow the skin to breathe...!

The affected parts of the body can be "streamed" once to five times a day with the AIRNERGY Stream for 5 to 20 minutes. For this purpose, the device can simply be placed on the affected body parts or held.

The AIRNERGY Stream can be used for the following complaints:

- Joint problems such as hip, knee, shoulder, elbow, toe, foot, finger, wrist and spinal pain
- Abdominal pain, headaches and migraines
- Bruises, sprains and burns
- Irritations of ears and sinuses, tinnitus and other ear complaints, eye problems, tooth and jaw pain as well as skin problems, cough, menstrual problems etc.

			
Power consumption	255 mA	Power consumption	510 mA
Power supply	5V DC 1.5 A	Power supply	5V DC 1.5 A
USB Power supply unit	100-240V AC 50-60Hz Out: 0.35 A	USB Power supply unit	100-240V AC 50-60Hz Out: 0.35 A
USB spiral cable	12V /0.3 A	USB spiral cable	12V /0.6 A
Dimensions (mm)	82 x115 x 50 (W x H x D)	Dimensions (mm)	239 x 124 x 285 (W x H x D)
Weight incl. bag	0.80 kg	Weight incl. bag	1.20 kg
			
Service	Maintenance free	Service	Maintenance free
Warranty	2 years	Warranty	2 years
Shipment includes	Stream USB power supply set, Transport bag, User manual	Shipment includes	Stream USB power supply set, transport bag, User manual, HS Diffusor
Article-No. (EAN-No.) silver carbon	3011055 (0792971115882) 3013055 (0792971115899)	Article-No. / EAN-No.	3054055 (2345678910125)

AIRNERGY_Dat_sheet_overview_2024